


Anti-Tead1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa, Chip - P28347

Item Number: CM0014290



Introduction

TEAD1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA, and ChIP with reactivity in human, mouse, and rat. Suitable for studying Hippo signaling, transcription regulation, and cardiac development.

[Learn More](#)

Parameter	Value
Product Name	TEAD1 Rabbit pAb
Remarks/Alias	AA; REF1; TCF13; TEF-1; NTEF-1; TCF-13; TEAD-1; TEAD1
Species	Human
GeneID (human)	7003
GeneID	7003
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human TEAD1 (NP_068780.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	P28347
Protein Weight	48kDa
Shipping	Ice bag