


Anti-Tead3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q99594

Item Number: CM0028006



Introduction

TEAD3 Rabbit pAb is a rabbit polyclonal antibody for WB, IHC-P, IF/ICC, ELISA, recognizing human, mouse, rat TEAD3 (TEF-5) at 49kDa. It targets the Hippo signaling transcription factor for placental and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TEAD3 Rabbit pAb
Remarks/Alias	TEF5; TEAD5; TEF-5; DTEF-1; ETFR-1; TEAD-3; TEAD3
Species	Human
GeneID (human)	7005
Immunogen	Recombinant fusion protein containing amino acids 130-330 of human TEAD3 (NP_003205.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99594
Protein Weight	49kDa
Shipping	Ice bag