

Anti-Tekt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q969V4

Item Number: CM0015229



Introduction

High-quality rabbit polyclonal antibody targeting human TEK1 (tektin-1), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying ciliary and flagellar microtubule biology.

[Learn More](#)

Parameter	Value
Product Name	TEKT1 Rabbit pAb
Remarks/Alias	TEKT1; tektin-1
Species	Human
GeneID (Human)	83659
GeneID	83659
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TEK1 (NP_444515.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q969V4
Protein Weight	48kDa
Shipping	Ice bag