

Anti-Tfam Rabbit Monoclonal Antibody For Wb, Ihc-P, If-P, Elisa - Q00059

Item Number: CM0007914



Introduction

Rabbit monoclonal antibody against human TFAM, validated for WB, IHC-P, IF-P, and ELISA. Cross-reacts with mouse and rat. Ideal for mitochondrial transcription and mtDNA maintenance research.

[Learn More](#)

Parameter	Value
Product Name	TFAM Rabbit mAb
Remarks/Alias	TCF6; MTTF1; MTTFA; TCF6L1; TCF6L2; TCF6L3; MTDPS15; TFAM
Species	Human
GeneID	7019
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human TFAM (NP_003192.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00059
Protein Weight	29kDa
Shipping	Ice bag