


Anti-Tfe3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P19532

Item Number: CM0003662



Introduction

Rabbit monoclonal antibody against human TFE3, validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Recombinant immunogen ensures high specificity for lysosomal biogenesis and immune response research.

[Learn More](#)

Parameter	Value
Product Name	TFE3 Rabbit mAb
Remarks/Alias	TFEA; RCCP2; RCCX1; MRXSPF; bHLHe33; TFE3
Species	Human
GeneID (Human)	7030
GeneID	7030
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 191-440 of human TFE3 (P19532).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19532
Protein Weight	62kDa
Shipping	Ice bag