

Anti-Tfr2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Up52

Item Number: CM0027997

Introduction

Rabbit polyclonal antibody against human TFR2 (Transferrin receptor protein 2). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse. For iron metabolism and hepatocyte research.

[Learn More](#)

Parameter	Value
Product Name	TFR2 Rabbit pAb
Remarks/Alias	HFE3; TFRC2; TFR2
Species	Human
Gene ID	7036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-330 of human TFR2 (NP_003218.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UP52
Protein Weight	89kDa
Shipping	Ice bag