

Anti-Tgs1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rs0

Item Number: CM0030502



Introduction

Anti-TGS1 rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human and mouse TGS1, a trimethylguanosine synthase involved in snRNA cap hypermethylation and Cajal body formation.

[Learn More](#)

Parameter	Value
Product Name	TGS1 Rabbit pAb
Remarks/Alias	PIMT; PIPMT; NCOA6IP; TGS1
Species	Human
GeneID (human)	96764
GeneID	96764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 714-853 of human TGS1 (NP_079107.6)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96RS0
Protein Weight	97kDa
Shipping	Ice bag