

Anti-Themis2 Polyclonal Antibody For Wb, Elisa - Q5Tej8

Item Number: CM0013995

Introduction

THEMIS2 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Ideal for macrophage inflammatory response and B cell activation research.

[Learn More](#)

Specification	Details
Product Name	THEMIS2 Rabbit pAb
Remarks/Alias	ICB1; ICB-1; C1orf38; THEMIS2
Species	Human
GeneID (Human)	9473
GeneID	9473
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human THEMIS2 (NP_004839.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5TEJ8
Protein Weight	72kDa
Shipping	Ice bag