


Anti-Thpo Polyclonal Antibody For Wb, Elisa - P40225

Item Number: CM0019673



Introduction

THPO Rabbit pAb is a polyclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat THPO. Detects a ~38 kDa protein; ideal for studying thrombopoiesis, platelet regulation, and megakaryocyte development.

[Learn More](#)

Parameter	Value
Product Name	THPO Rabbit pAb
Remarks/Alias	ML; TPO; MGDF; MKCSF; MPLLG; THCYT1; THPO
Species	Human
GeneID (Human)	7066
GeneID	7066
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-349 of human THPO (NP_001171068.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40225
Protein Weight	38kDa
Shipping	Ice bag