


Anti-Thymidine Kinase 1 (Tk1) Rabbit Monoclonal Antibody For Wb And Elisa - P04183

Item Number: CM0002438



Introduction

Rabbit monoclonal antibody targeting human Thymidine Kinase 1 (TK1), validated for WB and ELISA. Useful for studying nucleotide metabolism, cell cycle, and DNA synthesis. Recognizes a cytosolic enzyme with a predicted molecular weight of 25 kDa.

[Learn More](#)

Parameter	Value
Product Name	Thymidine Kinase 1 Rabbit mAb
Remarks/Alias	TK2; Thymidine Kinase 1
Species	Human
Gene ID (Human)	7083
Gene ID	7083
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Thymidine Kinase 1 (P04183).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04183
Protein Weight	25kDa
Shipping	Ice bag