


Anti-Tify10A Polyclonal Antibody For Wb, Elisa - Q9Lma8

Item Number: CM0019429



Introduction

TIFY10A Rabbit pAb (anti-JAZ1) is a polyclonal antibody for Western blotting and ELISA in *Arabidopsis thaliana*. Targets jasmonate signaling repressor; nucleus-localized. Ideal for plant defense research.

[Learn More](#)

Parameter	Value
Product Name	TIFY10A Rabbit pAb
Remarks/Alias	AtJAZ1; jasmonate-zim-domain protein 1; T29M8.5; T29M8_5; TIFY10A
Species	<i>Arabidopsis thaliana</i>
GeneID	838501
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 154-253 of <i>Arabidopsis thaliana</i> TIFY10A (NP_564075.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Arabidopsis thaliana</i>
SWISS	Q9LMA8
Protein Weight	28kDa
Shipping	Ice bag