

Anti-Tigd1/Eeyore Rabbit Polyclonal Antibody For Wb And Elisa - Q96Mw7

Item Number: CM0023003

Introduction

TIGD1 rabbit polyclonal antibody validated for WB and ELISA, reactive with human, mouse, rat. Targets TIGD1/Eeyore protein, a nuclear DNA-binding protein. Ideal for research on transposable element-derived proteins.

[Learn More](#)

Parameter	Value
Product Name	TIGD1 Rabbit pAb
Remarks/Alias	EEYORE; TIGD1
Species	Human
Gene ID (Human)	200765
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TIGD1 (NP_663748.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96MW7
Protein Weight	67kDa
Shipping	Ice bag