


Anti-Tigit Rabbit Polyclonal Antibody For Wb, Elisa - Q495A1

Item Number: CM0025123



Introduction

Validated for WB and ELISA, this anti-TIGIT rabbit polyclonal antibody recognizes human TIGIT and cross-reacts with mouse. Ideal for immune checkpoint research and T-cell immunoregulation studies.

[Learn More](#)

Parameter	Value
Product Name	TIGIT Rabbit pAb
Remarks/Alias	VSIG9; VSTM3; WUCAM
Species	Human
GeneID (Human)	201633
Immunogen	Recombinant fusion protein containing a.a. 22-141 of human TIGIT (NP_776160.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q495A1
Protein Weight	26kDa
Shipping	Ice bag