

Anti-Timm10B Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5J6

Item Number: CM0028523



Introduction

Rabbit polyclonal antibody targeting TIMM10B (FXC1/Tim9b), validated for WB, IHC-P, ELISA. Reacts with human and mouse. Ideal for studying mitochondrial protein import and TIM22 complex function.

[Learn More](#)

Parameter	Value
Product Name	TIMM10B Rabbit pAb
Remarks/Alias	FXC1; Tim9b; TIM10B; TIMM10B
Species	Human
GeneID (Human)	26515
GeneID	26515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-103 of human TIMM10B (NP_036324.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y5J6
Protein Weight	12kDa
Shipping	Ice bag