

Anti-Timm17A Polyclonal Antibody For Wb, If/Icc, Elisa - Q99595

Item Number: CM0030573

Introduction

Rabbit polyclonal antibody targeting mitochondrial import inner membrane translocase subunit Tim17-A. Validated for Western blot, IF/ICC, and ELISA in human and mouse samples. Swiss-Prot: Q99595.

[Learn More](#)

Parameter	Value
Product Name	TIMM17A Rabbit pAb
Remarks/Alias	TIM17; TIM17A; TIMM17A
Species	Human
Gene ID (Human)	10440
Gene ID	10440
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-171 of human TIMM17A (NP_006326.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99595
Protein Weight	18kDa
Shipping	Ice bag