

Anti-Timm22 Monoclonal Antibody For Wb And Elisa - Q9Y584

Item Number: CM0004677



Introduction

Rabbit monoclonal antibody against TIMM22 (Tim22), the voltage-gated channel subunit of the mitochondrial TIM22 translocase. Confirmed for WB and ELISA in human, mouse, and rat. Essential for multi-pass inner membrane protein import and mitochondrial disease research.

[Learn More](#)

Parameter	Value
Product Name	TIMM22 Rabbit mAb
Remarks/Alias	TEX4; TIM22; COXPD43; TIMM22
Species	Human
GeneID (Human)	29928
GeneID	29928
Immunogen	Synthetic peptide corresponding to amino acids 20-120 of human TIMM22
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y584
Protein Weight	20kDa
Shipping	Ice bag