

Anti-Timm29 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bsf4

Item Number: CM0022087



Introduction

High-quality rabbit polyclonal antibody targeting human TIMM29 (TIM29, C19orf52) for Western blot and ELISA applications. Validated in human and mouse samples, suitable for mitochondrial import studies.

[Learn More](#)

Parameter	Value
Product Name	TIMM29 Rabbit pAb
Remarks/Alias	TIM29; C19orf52; TIMM29
Species	Human
Gene ID (Human)	90580
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 78-260 of human TIMM29 (NP_612367.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BSF4
Protein Weight	29kDa
Shipping	Ice bag