

Anti-Tjp2 Polyclonal Antibody For Wb, Ip, Elisa - Q9Udy2

Item Number: CM0018549



Introduction

Rabbit polyclonal antibody against human TJP2 (ZO-2), validated for WB, IP, and ELISA. Cross-reacts with human and mouse. Targets tight junction protein ZO-2 involved in cell adhesion and signaling.

[Learn More](#)

Parameter	Value
Product Name	TJP2 Rabbit pAb
Remarks/Alias	ZO2; X104; FHCA1; PFIC4; DFNA51; DUP9q21.11; C9DUPq21.11; TJP2
Species	Human
GeneID (Human)	9414
GeneID	9414
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 951-1190 of human TJP2 (NP_004808.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UDY2
Protein Weight	134kDa
Shipping	Ice bag