

## Anti-Tle1 Polyclonal Antibody For Wb And Elisa - Q04724

Item Number: CM0014595



### Introduction

TLE1 Rabbit Polyclonal Antibody validated for WB and ELISA, cross-reacts with mouse. Ideal for studying transcriptional corepression in Wnt signaling and NF-κB pathways. Supplied by CamelBio.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | TLE1 Rabbit pAb                                                                                                    |
| Remarks/Alias    | ESG; ESG1; GRG1; TLE-1; TLE1                                                                                       |
| Species          | Human                                                                                                              |
| GeneID (human)   | 7088                                                                                                               |
| GeneID           | 7088                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 256-355 of human TLE1 (NP_005068.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Mouse                                                                                                              |
| SWISS            | Q04724                                                                                                             |
| Protein Weight   | 83kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |