

Anti-Tle2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q04725

Item Number: CM0004880

Introduction

Recombinant rabbit monoclonal antibody against human TLE2 (Transducin-like enhancer protein 2). Validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat samples. Suitable for Wnt signaling research.

[Learn More](#)

Parameter	Value
Product Name	TLE2 Rabbit mAb
Remarks/Alias	ESG; ESG2; GRG2; TLE2
Species	Human
Gene ID	7089
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-138 of human TLE2 (NP_003251.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04725
Protein Weight	80KDa
Shipping	Ice bag