

Anti-Tmed2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q15363

Item Number: CM0022154



Introduction

TMED2 rabbit polyclonal antibody for Western blot, immunofluorescence, and ELISA. Reacts with human, mouse, and rat TMED2, a key regulator of vesicular trafficking and GPCR exocytosis.

[Learn More](#)

Parameter	Value
Product Name	TMED2 Rabbit pAb
Remarks/Alias	p24; P24A; RNP24; p24b1; p24beta1; TMED2
Species	Human
Gene ID (Human)	10959
Gene ID	10959
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-170 of human TMED2 (NP_006806.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15363
Protein Weight	23kDa
Shipping	Ice bag