


Anti-Tnfr2/Tnfrsf1B Rabbit Polyclonal Antibody For Wb, Elisa - P20333

Item Number: CM0014301



Introduction

High-affinity anti-TNFR2/TNFRSF1B rabbit polyclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. Suitable for apoptosis and metabolic signaling research.

[Learn More](#)

Parameter	Value
Product Name	TNFR2/TNFRSF1B Rabbit pAb
Remarks/Alias	p75; TBPII; TNFBR; TNFR2; CD120b; TNFR1B; TNFR80; TNF-R75; p75TNFR; TNF-R-II; TNFR2/TNFRSF1B
Species	Human
Gene ID (Human)	7133
Gene ID	7133
Immunogen	Synthetic peptide corresponding to amino acids 295-377 of human TNFR2/TNFRSF1B (NP_001057.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P20333
Protein Weight	48kDa
Shipping	Ice bag