

Anti-Tnfrsf6B Polyclonal Antibody For Wb, Elisa - O95407

Item Number: CM0014587



Introduction

High-quality rabbit polyclonal antibody against human TNFRSF6B (Decoy receptor 3, DcR3). Validated for WB and ELISA, reacts with human and rat. Ideal for apoptosis and immune-signaling research.

[Learn More](#)

Parameter	Value
Product Name	TNFRSF6B Rabbit pAb
Remarks/Alias	M68; TR6; DCR3; M68E; DJ583P15.1.1; TNFRSF6B
Species	Human
GeneID (Human)	8771
GeneID	8771
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-300 of human TNFRSF6B (NP_003814.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O95407
Protein Weight	33kDa
Shipping	Ice bag