

Anti-Tomm22 Polyclonal Antibody For Wb And Elisa - Q9Ns69

Item Number: CM0015542



Introduction

Rabbit polyclonal antibody targeting human TOMM22, a mitochondrial import receptor subunit, validated for WB and ELISA. Ideal for mitochondrial import, mitophagy, and preprotein translocation research.

[Learn More](#)

Parameter	Value
Product Name	TOMM22 Rabbit pAb
Remarks/Alias	1C9-2; TOM22; MST065; MSTP065; TOMM22
Species	Human
GeneID (Human)	56993
GeneID	56993
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-83 of human TOMM22 (NP_064628.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NS69
Protein Weight	16kDa
Shipping	Ice bag