

Anti-Topors Rabbit Polyclonal Antibody For Wb And Elisa - Q80Z37

Item Number: CM0019144

Introduction

Rabbit polyclonal antibody against Topors (LUN/p53BP3), an E3 ubiquitin and SUMO1 ligase involved in tumor suppression. Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for studies on p53 regulation and DNA damage response.

[Learn More](#)

Parameter	Value
Product Name	Topors Rabbit pAb
Remarks/Alias	LUN; TP53BP3; p53BP3/LUN; Topors
Species	Human
GeneID (Human)	106021
GeneID	106021
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of mouse Topors (NP_598858.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q80Z37
Protein Weight	117kDa
Shipping	Ice bag