

Anti-Tph1 Polyclonal Antibody For Wb, Elisa - P17752

Item Number: CM0022031



Introduction

TPH1 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, recognizing human, mouse, and rat TPH1. Ideal for serotonin biosynthesis research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	TPH1 Rabbit pAb
Remarks/Alias	TPRH; TRPH; TPH1
Species	Human
GeneID (Human)	7166
GeneID	7166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 345-444 of human TPH1 (NP_004170.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17752
Protein Weight	51kDa
Shipping	Ice bag