

Anti-Tph1 Rabbit Monoclonal Antibody For Ihc-P, Elisa - P17752

Item Number: CM0010978



Introduction

High-quality rabbit monoclonal antibody against human TPH1 (tryptophan hydroxylase 1). Validated for IF, IHC-P, and ELISA with cross-reactivity to mouse and rat. Ideal for serotonin biosynthesis studies.

[Learn More](#)

Parameter	Value
Product Name	TPH-1 Rabbit mAb
Remarks/Alias	TPRH; TRPH
Species	Human
GeneID (human)	7166
GeneID	7166
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 102-402 of human TPH1 (NP_004170.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-F, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17752
Protein Weight	51kDa/53kDa
Shipping	Ice bag