


Anti-Traf6 Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Y4K3

Item Number: CM0007450



Introduction

Rabbit monoclonal antibody against human TRAF6, validated for WB, IHC-P, IP, and ELISA. Reacts with human, mouse, and rat. Recognizes TRAF6, an E3 ubiquitin ligase crucial for NF- κ B and JNK signaling in immunity and osteoclast differentiation. Suitable for research and assay development.

[Learn More](#)

Parameter	Specification
Product Name	TRAF6 Rabbit mAb
Remarks/Alias	RNF85; MGC:3310; TRAF6
Species	Human
GeneID (Human)	7189
GeneID	7189
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4K3
Protein Weight	60 kDa
Shipping	Ice bag