


Anti-Traf6 Rabbit Polyclonal Antibody For If/Icc And Elisa - Q9Y4K3

Item Number: CM0016826



Introduction

Rabbit polyclonal antibody against human TRAF6 (SWISS Q9Y4K3), suitable for IF/ICC and ELISA. Cross-reacts with mouse and rat. Validated for research use in immune signaling and osteoclast differentiation.

[Learn More](#)

Parameter	Value
Product Name	TRAF6 Rabbit pAb
Remarks/Alias	RNF85; MGC:3310; TRAF6
Species	Human
Gene ID	7189
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 50-150 of human TRAF6 (NP_004611.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4K3
Protein Weight	60kDa
Shipping	Ice bag