


Anti-Transcription Factor Ap-2-Alpha (Tfap2A) Monoclonal Antibody For Wb, Elisa - P05549

Item Number: CM0007491



Introduction

Rabbit monoclonal antibody raised against a synthetic peptide within amino acids 50-150 of human Transcription factor AP-2-alpha (P05549). Validated for WB, ELISA applications; cross-reacts with human and mouse samples. Ideal for transcription regulation and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Activator protein 2 (AP-2/TFAP2A) Rabbit mAb
Remarks/Alias	AP-2; BOFS; AP2TF; TFAP2; AP-2alpha; Activator protein 2 (AP-2/TFAP2A)
Species	Human
Gene ID (Human)	7020
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Activator protein 2 (AP-2/TFAP2A) (P05549).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05549
Protein Weight	48kDa
Shipping	Ice bag