

Anti-Trappc10 Polyclonal Antibody For Wb, Elisa - P48553

Item Number: CM0030642

Introduction

Rabbit polyclonal antibody targeting human TRAPPC10, a TRAPP II complex subunit, validated for WB and ELISA. Cross-reacts with mouse and rat, suitable for Golgi trafficking and neurological disease research.

[Learn More](#)

Parameter	Value
Product Name	TRAPPC10 Rabbit pAb
Remarks/Alias	EHOC1; GT334; TMEM1; TRS30; EHOC-1; TRS130; NEDMISS; TRAPPC10
Species	Human
GeneID (Human)	7109
GeneID	7109
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human TRAPPC10 (NP_003265.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P48553
Protein Weight	142kDa
Shipping	Ice bag