


## Anti-Trem2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nzc2

Item Number: CM0014586



### Introduction

High-quality rabbit polyclonal antibody targeting human TREM2 (Q9NZC2). Validated for WB, IF/ICC, ELISA, cross-reacts with mouse and rat. Ideal for microglial biology and Alzheimer's research.

[Learn More](#)

| Parameter        | Value                                                                                                                                  |
|------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TREM2 Rabbit pAb                                                                                                                       |
| Remarks/Alias    | PLOSL2; TREM-2; Trem2a; Trem2b; Trem2c; TREM2                                                                                          |
| Species          | Human                                                                                                                                  |
| Gene ID (Human)  | 54209                                                                                                                                  |
| Gene ID          | 54209                                                                                                                                  |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-161 of human TREM2 (NP_061838.1). |
| Source           | Rabbit                                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                                  |
| Application      | WB, IF/ICC, ELISA                                                                                                                      |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                      |
| SWISS            | Q9NZC2                                                                                                                                 |
| Protein Weight   | 25kDa                                                                                                                                  |
| Shipping         | Ice bag                                                                                                                                |