

Anti-Trif/Ticam1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8IUC6

Item Number: CM0016005



Introduction

Recombinant-based rabbit polyclonal antibody targeting human TRIF/TICAM1 (UniProt Q8IUC6). Validated for WB, ELISA, and cross-reacts with mouse and rat. Ideal for TLR signaling and innate immunity research.

[Learn More](#)

Parameter	Information
Product Name	TRIF/TICAM1 Rabbit pAb
Remarks/Alias	TRIF; IIAE6; MyD88-3; PRVTIRB; TICAM-1; TRIF/TICAM1
Species	Human
Gene ID (Human)	148022
Gene ID	148022
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-200 of human TRIF/TICAM1 (NP_891549.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8IUC6
Protein Weight	76kDa
Shipping	Ice bag