

Anti-Trim23 Polyclonal Antibody For Wb, If/Icc, Elisa - P36406

Item Number: CM0030403



Introduction

Rabbit polyclonal antibody against human TRIM23 (ARD1/ARFD1/RNF46), validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studying autophagy, ubiquitination, and viral infection pathways.

[Learn More](#)

Parameter	Value
Product Name	TRIM23 Rabbit pAb
Remarks/Alias	ARD1, ARFD1, RNF46, TRIM23
Species	Human
GeneID (Human)	373
GeneID	373
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human TRIM23 (NP_001647.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36406
Protein Weight	64kDa
Shipping	Ice bag