


Anti-Trim9 Rabbit Polyclonal Antibody For Wb, Elisa - Q9C026

Item Number: CM0013753



Introduction

TRIM9 Rabbit Polyclonal Antibody validated for Western Blot and ELISA, showing reactivity with human, mouse, and rat TRIM9 (E3 ubiquitin ligase). Ideal for neuroscience and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	TRIM9 Rabbit pAb
Remarks/Alias	RNF91, SPRING, TRIM9
Species	Human
GeneID (human)	114088
GeneID	114088
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human TRIM9 (NP_055978.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9C026
Protein Weight	79kDa
Shipping	Ice bag