

Anti-Trmt1 Polyclonal Antibody For Wb, Elisa - Q9Nxx9

Item Number: CM0026863



Introduction

Rabbit polyclonal antibody against human tRNA methyltransferase TRMT1 (TRM1; MRT68). Validated for WB and ELISA; reacts with human and mouse. Suitable for tRNA modification and redox biology research.

[Learn More](#)

Parameter	Value
Product Name	TRMT1 Rabbit pAb
Remarks/Alias	TRM1; MRT68; TRMT1
Species	Human
Gene ID (Human)	55621
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NXH9
Protein Weight	72kDa
Shipping	Ice bag