


Anti-Trmt2A/Htf9C Rabbit Monoclonal Antibody For Wb, Elisa - Q8Iz69

Item Number: CM0003676



Introduction

Rabbit monoclonal antibody targeting human TRMT2A/HTF9C, a tRNA/mRNA methyltransferase. Validated in WB and ELISA. Detects a band at ~69 kDa. Suitable for RNA modification studies.

[Learn More](#)

Parameter	Value
Product Name	TRMT2A/HTF9C Rabbit mAb
Remarks/Alias	HTF9C; TRMT2A/HTF9C
Species	Human
GeneID	27037
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TRMT2A/HTF9C (Q8Iz69)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8Iz69
Protein Weight	69kDa
Shipping	Ice bag