

Anti-Tropomodulin-2 (Tmod2) Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q9Nzr1

Item Number: CM0002771

Introduction

Rabbit monoclonal antibody targeting human TMOD2 (tropomodulin-2), a neuronal actin-capping protein. Validated for WB, IF-P, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for cytoskeleton and neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	TMOD2 Rabbit mAb
Remarks/Alias	NTMOD; N-TMOD; TMOD2
Species	Human
Gene ID	29767
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human TMOD2 (Q9NZR1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NZR1
Protein Weight	40kDa
Shipping	Ice bag