

Anti-Troponin I2 (Tnni2) Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P48788

Item Number: CM0008829



Introduction

CamelBio offers rabbit monoclonal Troponin I2 (TNNI2) antibody validated for WB, IF-P, ELISA. Detects human, mouse, and rat fast skeletal muscle troponin I. 21 kDa. Ideal for muscle contraction and myopathy research.

[Learn More](#)

Parameter	Value
Product Name	Troponin I2 (TNNI2) Rabbit mAb
Remarks/Alias	DA2B, FSSV, DA2B1, fsTnl, AMCD2B, Troponin I2 (TNNI2)
Species	Human
GeneID (Human)	7136
GeneID	7136
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 83-182 of human Troponin I2 (TNNI2) (P48788)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P48788
Protein Weight	21 kDa
Shipping	Ice bag