


Anti-Trpa1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - O75762

Item Number: CM0028021



Introduction

TRPA1 Rabbit pAb is a polyclonal antibody targeting human TRPA1, suitable for WB, IF-P, IHC-P, and ELISA. Cross-reacts with mouse and rat. Useful for pain and sensory transduction research.

[Learn More](#)

Parameter	Value
Product Name	TRPA1 Rabbit pAb
Remarks/Alias	FEPS; p120; FEPS1; ANKTM1; TRPA1
Species	Human
GeneID (Human)	8989
GeneID	8989
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human TRPA1 (NP_015628.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75762
Protein Weight	128 kDa
Shipping	Ice bag