


Anti-Trpm8 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z2W7

Item Number: CM0029724



Introduction

High-quality TRPM8 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human TRPM8, a cold-sensing ion channel implicated in thermosensation and prostate cancer. Ideal for research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	TRPM8 Rabbit pAb
Remarks/Alias	TRPP8; LTRPC6; trp-p8; LTrpC-6; TRPM8
Species	Human
GeneID (Human)	79054
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 980-1104 of human TRPM8 (NP_076985.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z2W7
Protein Weight	128kDa
Shipping	Ice bag