


Anti-Trpv5 Polyclonal Antibody (Rabbit) For Wb, Ihc-P, Elisa - Q9Nqa5

Item Number: CM0025872



Introduction

High-quality rabbit polyclonal antibody against human TRPV5 (CAT2/ECAC1). Validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for calcium transport and kidney research.

[Learn More](#)

Parameter	Value
Product Name	TRPV5 Rabbit pAb
Remarks/Alias	CAT2; ECAC1; OTRPC3; TRPV5
Species	Human
GeneID (Human)	56302
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 655-725 of human TRPV5 (NP_062815.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQA5
Protein Weight	83kDa
Shipping	Ice bag