

## Anti-Trβ1/Thrb Polyclonal Antibody For Wb, Elisa - P10828

Item Number: CM0020158



### Introduction

Rabbit polyclonal antibody targeting human thyroid hormone receptor beta (TRβ1, THRB) for Western blotting and ELISA. Cross-reacts with mouse and rat. Suitable for nuclear receptor research, transcription regulation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                             |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TRβ1/THRB Rabbit pAb                                                                                                              |
| Remarks/Alias    | TRb; GRTH; PRTH; THR1; ERBA2; NR1A2; THRB1; THRB2; TRbeta; THRbeta; TRbeta1; C-ERBA-2; THRbeta1; Thrbeta2; C-ERBA-BETA; TRβ1/THRB |
| Species          | Human                                                                                                                             |
| GeneID (Human)   | 7068                                                                                                                              |
| GeneID           | 7068                                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-263 of human TRβ1/THRB (NP_000452.2)              |
| Source           | Rabbit                                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                                             |
| Application      | WB, ELISA                                                                                                                         |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                 |
| SWISS            | P10828                                                                                                                            |
| Protein Weight   | 53kDa                                                                                                                             |
| Shipping         | Ice bag                                                                                                                           |