

Anti-Tryptophanyl-Trna Synthetase 1 Polyclonal Antibody For Wb, Elisa - P23381

Item Number: CM0029005



Introduction

Rabbit polyclonal antibody targeting human Tryptophanyl-tRNA synthetase 1 (TrpRS/WARS). Validated for WB and ELISA, cross-reacts with mouse and rat. Suitable for protein translation and angiogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	Tryptophanyl-tRNA synthetase 1 Rabbit pAb
Remarks/Alias	HMN9; WARS; IFI53; IFP53; GAMMA-2; Tryptophanyl-tRNA synthetase 1
Species	Human
GeneID (Human)	7453
GeneID	7453
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human Tryptophanyl-tRNA synthetase 1 (NP_776049.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23381
Protein Weight	53kDa
Shipping	Ice bag