

Anti-Tspan18 Polyclonal Antibody For Wb, Elisa - Q80Wr1

Item Number: CM0024013

Introduction

Rabbit polyclonal antibody against mouse Tspan18, suitable for Western blot and ELISA applications, validated for mouse and rat reactivity. Tspan18 regulates ORAI1 localization and Ca²⁺ signaling.

[Learn More](#)

Parameter	Value
Product Name	Tspan18 Rabbit pAb
Remarks/Alias	6720430015; 2610042G18Rik; Tspan18
Species	Mouse
GeneID	241556
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 149-248 of mouse Tspan18 (NP_899003.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q80WR1
Protein Weight	28kDa
Shipping	Ice bag