

Anti-Tspan8 Rabbit Polyclonal Antibody For Wb, Elisa - P19075

Item Number: CM0013717

Introduction

Anti-TSPAN8 rabbit polyclonal antibody validated for WB and ELISA. Cross-reacts with human, mouse, and rat. Target tetraspanin-8 involved in cell signaling, migration, and tumor progression. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	TSPAN8 Rabbit pAb
Remarks/Alias	CO-029; TM4SF3; TSPAN8
Species	Human
GeneID (Human)	7103
GeneID	7103
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-200 of human TSPAN8 (NP_004607.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P19075
Protein Weight	26kDa
Shipping	Ice bag