

Anti-Tspy3 Rabbit Polyclonal Antibody For Wb And Elisa - P0Cv98

Item Number: CM0019444

Introduction

TSPY3 rabbit polyclonal antibody validated for WB and ELISA. Detects human TSPY3, cross-reacts with rat. For sperm differentiation and proliferation research. - P0CV98

[Learn More](#)

Parameter	Value
Product Name	TSPY3 Rabbit pAb
Remarks/Alias	CT78; TSPY; TSPY1; TSPY10; TSPY3
Species	Human
GeneID (Human)	728137
GeneID	728137
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-308 of human TSPY3 (NP_001071165.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P0CV98
Protein Weight	35kDa
Shipping	Ice bag