


Anti-Uba6 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - A0Avt1

Item Number: CM0026263



Introduction

Validated rabbit polyclonal antibody against UBA6 (UBE1L2, E1-L2), a ubiquitin-activating enzyme. Suitable for WB, IF/ICC, IP, ELISA on human and mouse samples. Ideal for ubiquitin pathway and fertility research.

[Learn More](#)

Parameter	Value
Product Name	UBA6 Rabbit pAb
Remarks/Alias	E1-L2, MOP-4, UBE1L2, UBA6
Species	Human
GeneID	55236
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human UBA6 (NP_060697.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	A0AVT1
Protein Weight	118kDa
Shipping	Ice bag