

Anti-Ube2L6 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O14933

Item Number: CM0006633



Introduction

Anti-UBE2L6 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat UBE2L6 (UBCH8/RIG-B), an E2 enzyme in ubiquitin and ISG15 conjugation. Useful for ubiquitin-proteasome research.

[Learn More](#)

Parameter	Value
Product Name	UBE2L6 Rabbit mAb
Remarks/Alias	RIG-B; UBCH8; UBE2L6
Species	Human
Gene ID	9246
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 54-153 of human UBE2L6 (NP_004214.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14933
Protein Weight	18kDa
Shipping	Ice bag