

Anti-Ube2S Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q16763

Item Number: CM0026899



Introduction

High-specificity rabbit polyclonal antibody targeting human UBE2S, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying APC/C-mediated ubiquitination and mitotic progression.

[Learn More](#)

Parameter	Value
Product Name	UBE2S Rabbit pAb
Remarks/Alias	EPF5; E2EPF; E2-EPF; UBE2S
Species	Human
GeneID (Human)	27338
GeneID	27338
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human UBE2S (NP_055316.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16763
Protein Weight	24kDa
Shipping	Ice bag